



GE Light Speed 16

The LightSpeed 16 CT Scanner features an innovative solid state LightSpeed Matrix detector that enables the system to simultaneously acquire multiple adjacent slices of the anatomy. LightSpeed systems provide exceptionally clear image quality, short scan times and accurate non-invasive diagnosis for improved patient outcome. The underlying principles are the powerful 53 kW generator, the 6.3 mhu tube and the 500 millisecond scan time. In tandem with the high detection efficiency of the LightSpeed Detector, these components provide high signal output at very low patient dose.



- 36 Months \$ 10,730.30 The cost to maintain this GE Light Speed 16 is \$298.00 per day
- 60 Months \$ 8,266.26 The cost to maintain the GE Light Speed 16 is \$137.00 per day

(The cost to maintain the GE Light Speed doesn't include electric, Ct Tech, Insurance or rent)

**Have a question
or need to know more?**

Click here to contact us



*or Call Now at
732-262-3115*

NJ Warehouse
1111 Industrial Pkwy
Brick, NJ 08724
Tel: (732) 262-3115

NJ Headquarters
638 Mantoloking Road
Brick, NJ 08723
Tel: (732) 262-3115 Fax: (732) 262-3105

Florida Office
11176 Creek Haven Dr.
Riverview, FL 33569
Tel: (732) 262-3115

Email: Sales@nationwideimaging.com

www.NationwideImaging.com